

PHILIPS

CarePoint



Pull Cord

Pull cord is ideal for bathrooms and bedsides where help is frequently needed. Pull cords are available either for flush or surface mount.

UL2560
Certified

Features & Benefits

Magnetic Switch is ideal for bathrooms, where steam and spray vapors may exist

Wireless Supervised Transmitter

Available as Flush Mount or Surface Mount

Long-Life 3 Volt Lithium battery changeable on-site

Resetting 'Emergency' switch automatically acknowledges Alarm

Supervised device provides notification when device needs service or battery replacement

Specifications

Color	White with red slide switch
Cord	6ft. white nylon with white faub
Dimensions	Surface Mount: 3.125" W x 4.75" H x 2" D
	Flush Mount: 2.75" W x 4.5" H x 1" D
Battery	3V DC Lithium
Average Battery Life	5 years
Compliance/Regulatory	UL2560, FCC Part 15

Mounting Note: Flush Mount devices require a non-metallic single 18.0 cu. in. gang box

