



Passive Infared (PIR) Motion Sensor

Passive check-in device: when movement is detected, a signal is sent to the Central Monitoring Station documenting the activity.

Features & Benefits

Can be utilized for automated resident check-in or activity monitoring

A wall tamper triggers an alarm if the sensor is removed

3-minute lockout timer extends battery life

Wall-mounting bracket

Supervised device provides notification when device needs service or battery replacement

Specifications

Color	White
Dimensions	2.75" W x 4.38" H x 1.75" D
Battery (1)	3.6V DC Lithium
Average Battery Life	4-5 years
Installation	Wall mounting



www.philips.com