

PHILIPS

CarePoint



Resident Check-In

Active check-in device: when a resident touches it, a signal is sent to the Central Monitoring Station documenting the activity.

Features & Benefits

Sends resident check-in messages to the Central Monitoring Station

Reports check-ins to the SmartCare Software allowing report generation noting residents which have not yet checked in

Activation is simple, press anywhere on the face for resident check-in

Easy-to-clean membrane switch

Wall mounted for resident convenience

Supervised device provides notification when device needs service or battery replacement

Specifications

Color	White with red face
Dimensions	2.56" W x 4" H x .94" D
Batteries	3.6V DC Lithium
Average Battery Life	5 years

© 2015. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. (Royal Philips) or their respective owners.

www.philips.com



111 Lawrence Street
Framingham, MA 01702
Telephone: 800-816-4885

www.philipsseniorliving.com