

**PHILIPS**

CarePoint



# Roam Alert® Wandering Tag

Roam Alert® Wandering Tags provide protection to wandering-prone residents.

The Wandering Tag can be worn on the individual's wrist or ankle and is attached using a tear-proof band. The tag is waterproof and is made with a smooth surface that permits easy sterilization and prevents harboring of bacteria.

## Specifications

<b>Output Frequency</b>	433.92 MHz
<b>Input Frequency</b>	307 kHz
<b>Battery Life</b>	Up to three years
<b>Operating Temperature</b>	-10° C to 50° C (14° F to 122° F)
<b>Humidity</b>	0-95% RH @ 70° F (21° C), non-condensing
<b>Dimensions</b>	Approximately 1.1 x 1.0 x 0.5 in (2.79 x 2.54 x 1.27 cm)
<b>Weight</b>	Approximately 0.3 oz. (8.5 g)

## Ordering Information

<b>Description</b>	<b>Part Number</b>
RoamAlert Wrist Tags	1900547
Roam Alert Tag Rack	1900551

© 2015. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. (Royal Philips) or their respective owners.

[www.philips.com](http://www.philips.com)



111 Lawrence Street  
Framingham, MA 01702  
Telephone: 800-816-4885

[www.philipsseniorliving.com](http://www.philipsseniorliving.com)